

HTL tips	GILSON pipettes		Pipetman Classic single channel						Pipetman Neo single channel						Pipetman Neo multi		Pipetman Ultra single channel						Pipetman Ultra multi		
	Tip volume [µl]	Cat number	0,2-2	1-10	2-20	20-100	50-200	200-1000	0,2-2	1-10	2-20	20-100	20-200	200-1000	1-20	20-300	0,2-2	1-10	2-20	20-100	20-200	200-1000	1-20	20-300	
NON FILTER	0,5-10	30011	•	•					•	•							•	•							
		30012	•	•					•	•								•	•						
		30013	•	•					•	•								•	•						
		30014	•	•					•	•								•	•						
		30015	•	•					•	•								•	•						
	1-200	30201 / 30221			•	•	•				•	•	•						•	•	•				
		30202 / 30222			•	•	•				•	•	•						•	•	•				
		30203 / 30223			•	•	•				•	•	•						•	•	•				
		30224			•	•	•				•	•	•						•	•	•				
		30225			•	•	•				•	•	•						•	•	•				
	1-300	30301				•	•					•	•			•				•	•				•
		30302				•	•					•	•			•				•	•				•
		30303				•	•					•	•			•				•	•				•
	100-1000	31001						•						•									•		
		31002						•						•									•		
		31003						•						•									•		
	100-1200	31201						•						•									•		
		31202						•						•									•		
	FILTER	0,5-10	30016	•	•					•	•							•	•						
		0,5-20	30026			•					•								•						
0,5-50		30056																							
1-100		30106				•					•								•						
1-200		30206					•					•			•					•				•	
1-300		30306																							
100-1000		31006						•					•									•			
100-1100		31206						•					•									•			